


# Experiment P007: Acceleration due to Gravity (Free Fall Adapter)



## EQUIPMENT NEEDED

*Science Workshop™ Interface*  
Clamp, right angle  
Base and support rod  
Free fall adapter  
Balls, 13 mm and 19 mm  
Meter stick (or metric tape measure)

## PURPOSE

The purpose of this laboratory activity is to measure the acceleration of a falling object assuming that the only force acting on the object is the gravitational force.

## THEORY

One equation describing the motion of a body starting from rest and undergoing constant acceleration can be expressed as:

$$d = \frac{1}{2}at^2$$

where **d** is the distance the object has traveled from its starting point, **a=g** is the acceleration of the object, and **t** is the time elapsed since the motion began.

Therefore, the acceleration is equal to the distance doubled, divided by the time squared, or:

$$a=2d/t^2$$

In order to measure the acceleration caused by the gravitational force, several questions must be answered:

Is the acceleration constant? If it is, then the distance the object falls will be proportional to the square of the elapsed time, as in the above equation.

If the acceleration is constant, what is the value of the acceleration? Is it the same for all objects or does it vary with mass, size, or with some other property of the object? If it is not constant, how does it vary with time?

## **PROCEDURE**

For this laboratory activity, the free fall adapter will measure the time of fall for steel balls of different sizes and mass. Using a meter stick, you will directly measure the distance the ball falls. The *DataStudio* program will then plot distance doubled vs. the time squared. Then the slope of the (best fit line) for this graph will then give the acceleration due to gravity, *g*.

### **PART I: Computer Setup**

1. Connect the *Science Workshop* interface to the computer, turn on the interface, and turn on the computer.
2. Connect the free fall adapter's stereo phone plug into Digital Channel 1.
3. Open the *DataStudio* document titled as shown:

Windows: P007\_ADAP.DS

The document will open with a Graph display of Distance "(m)" versus "Time (sec)", and a Table display of "Free Fall Data", the time of fall.

The Experimental Setup window has been resized. If you want to expand the Experiment Setup window to its original size, click on the "Zoom" box in the upper right hand corner of the window. (Note: To bring a display to the top, click on its window or select the name of the display from the list at the end of the Display menu.)

## **PART II: Sensor Calibration and Equipment Setup**

You do not need to calibrate the free fall adapter.

1. Attach the right angle clamp to a vertical support rod. Mount the free fall adapter's release mechanism horizontally in the clamp.
2. Place the free fall adapter's timing pad on the floor directly below the release mechanism.
3. Place the larger ball (19 mm) in the release mechanism. Press the spring-loaded rod inward to hold the ball in the mechanism, and tighten the thumbscrew to hold the rod in place.
4. Adjust the position of the release mechanism so the bottom of the ball is 1.50 meters above the timing pad.
5. Practice dropping the ball a few times before taking data.

## **PART IIIA: Data Recording – 19 mm Ball**

1. Measure the distance from the bottom of the ball to the top of the timing pad. Record this distance in the Free Fall Data table.
2. Click the “Start” button to begin data recording.
3. Loosen the thumbscrew to release the ball.
4. After the ball hits the timing pad, click “Keep”, the Manual Sampling window will open, type in the distance the ball fell. Click “OK” to record your distance and time.

The distance you typed in appears in the Free Fall Data Table and the default value changes.

5. Once again place the ball in the release mechanism. Adjust the position of the release mechanism so the bottom of the ball is 1.25 meters above the timing pad. Measure the distance between the ball and the pad and record it in the Free Fall Data Table.
6. When you are ready to record the next drop, loosen the thumbscrew to release the ball.
7. After the ball hits the pad, click “Keep”, type in the new distance (“1.25”) in the Manual Sampling window. Click “OK” to record the distance and time in the data table below. Your new distance appears in the Free Fall Data table and the default value changes.

8. Replace the ball in the mechanism. Move the mechanism so the bottom of the ball is 1.00 meters from the top of the timing pad. Measure and record the distance.
9. Loosen the thumbscrew to release the ball.
10. After the ball hits the timing pad, type in the distance the ball fell in the Manual Sampling window. Click “OK” to record your distance and time in the table below.
11. Repeat the process for heights of 0.75 and 0.50 meters.
12. Click the “Stop” button to end data recording.
13. The Manual Sampling window will close. Run #1 will appear in the Data list under Time of Fall, Ch 1(s) and under Distance (m).
14. Click the Graph to make it active. Select “Linear” from the “Fit” menu.

The slope of the best-fit line is the acceleration due to gravity,  $g$ , for the most recent run of data.

**DATA TABLE: 19 mm Ball**

Trial	1	2	3	4	5
Distance (m)					
Time (sec)					

**Acceleration due to gravity (slope) = \_\_\_\_\_ m/sec/sec**

### PART IIIB: Data Recording – 13 mm Ball

1. Place the smaller ball (13 mm) in the release mechanism. Adjust the position of the release mechanism so the bottom of the ball is 1.50 meters above the timing pad. Measure the distance from the bottom of the 13 mm ball to the top of the timing pad. Record this distance in the Free Fall Data table.
2. Click the “START” button to begin data recording.
3. Loosen the thumbscrew to release the ball.
4. After the ball hits the timing pad, type in the distance the ball fell in the Manual Sampling window. Click “OK” to record your distance and time.
5. Repeat the data recording process with the 13 mm ball for distances of 1.25, 1.00, 0.75, and 0.50 meters as in Part IIIA.
6. When all the distances have been measured, click the “Stop” button to end data recording.
7. The Manual Sampling window will close. Run #2 will appear in the Data list under Time of Fall, Ch 1(s) and under Distance (m).
8. Click the Graph to make it active. Select “Linear” from the “Fit” menu.

The slope of the best-fit line is the acceleration due to gravity,  $g$ , for the most recent run of data.

#### DATA TABLE: 13 mm Ball

Trial	1	2	3	4	5
Distance (m)					
Time (sec)					

Acceleration due to gravity (slope) = \_\_\_\_\_ m/sec/sec

## QUESTIONS

1. What is the percent difference between the measured value for the acceleration due to gravity and the accepted value (9.8 m/sec/sec)?

	<b>% Difference</b>
<b>19 mm ball</b>	
<b>13 mm ball</b>	

2. What are possible reasons for the difference in the accepted value of the acceleration of a free falling object,  $g = 9.8 \text{ (m/s}^2\text{)}$ .
3. Is the experimental value of acceleration due to gravity the same for both the 19 mm ball and the 13 mm ball?
4. What are possible reasons for this difference, if any?